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Warranty Card

Congratulations on your Sensica purchase!

We have every faith in the quality of each Sensica device - we know they meet the highest standards, because we subject them to rigorous quality and performance testing. In the unlikely event that you do find a defect in a Sensica device, it may be covered by this warranty.

Two Year Limited Warranty

This warranty covers only Sensica devices. Other Sensica products, such as pads or creams, are defined as consumables and therefore excluded from warranty coverage.This limited warranty covers manufacturer defect related problems in your Sensica device, if you purchased it directly from www.sensica.com or from an authorized partner. Warranty approvals will only be granted where an original receipt is provided. This warranty is limited to devices proven to be defective. In order to obtain warranty service on your Sensica device, please go to our website or call our customer service center for assistance (you will find the number at www.sensica.com). Sensica requires proof of purchase. Please provide an invoice or purchase receipt showing the date of purchase of the device.

Enter your serial number on our website
to activate your warranty.

www.sensica.com

Model: ST100

Double insulation

Consult Instructions for use

Waste Electronic and Electrical Equipment

Please read all warnings and safety precautions prior to use

Protect from moisture

This device must not be used near to water or in water itself (e.g., in a wash basin, in the shower or bathtub) – risk of electric shock!

Our products and packaging are recyclable, do not throw them away!

CE

UK CA

Use the device only with the original power supply. For your safety, never use DC adapters which are not approved by the manufacturer.

English

MD 04565 A0 / RD 11012 H0

Device specification parameters

The RF handheld device is powered by bipolar RF at 10±0.5Watt

Frequency
1±0.05MHz

Selectable output power level
3 levels

Mechanical massage
Implemented by moving RF electrodes.

Speed
50-100 RPM

Temperature sensors
RF heating is automatically deactivated when skin temperature reaches 40°C/104°F. The heating is reactivated when temperature decreases.

Treatment area
~2.2cm² / 0.34in²

External power supply
Input 100-240VAC, 50-60Hz,
Output 9VDC, 1.5A

Device size and weight

140g, 173mm x 56.8mm x 53.1mm
4.93oz, 6.8in x 2.23in x 2.09in

Storage Instructions

Storage and transport conditions	
Temperature	-25°C to 70°C (-13°F to 158°F)
Humidity	10% to 90% RH at 55°C/131°F
Pressure	Non-condensing
Operational conditions	
Temperature	5°C to 40°C (41°F to 104°F)
Humidity	20% to 80% RH
Max operating altitude	3,000 m.

If you cannot feel heat:

- Verify full contact with the skin (indicated by a green light).
- Verify that you applied enough Base Gel.
- Remember that the device reacts to your skin and may temporarily reduce heat for safety reasons.

See status indicator table in the **User Interface Description** section.

In case of dirt or debris on the treatment fingers:

- you should clean them in accordance with the **Cleaning & Storage** section. If you can't get the area clean, contact our customer service.

In case of a general error:

- reset the device by unplugging the power cord from the electrical outlet on the wall, wait 10 seconds and then reconnect it.

If you still have problems with your device, please contact our customer service.

Customer Service

If your device has broken/ stopped working or if you need any further assistance/ information, please contact us. You may also contact us to order additional Base Gel or you can purchase it online at sensica.com.

For customer service, please see contact details on inner back cover.

Warning: avoid contact with the eyes, eye cavities, nostrils and oral cavity, as treatment in these areas may cause adverse effects.

Warning: avoid treatment in the lower third of the neck since this area is more sensitive and avoid treating the thyroid area (marked ⚠).

- Response to the treatment differs between individuals and depends on environmental and biological factors, heredity and aging.

Device Maintenance

Cleaning & Storage

- Make sure to turn off and unplug the device before cleaning.
- It is recommended to clean the device, especially the treatment fingers, temperature sensors, and silicone cover, of residual Base Gel with a soft damp cloth or wipes.

Do not rinse the device under running water, as it may cause device malfunction.

- Keep the device dry. Never place, store or use the device near: bathtubs, showers, basins or other vessels containing water as this may cause severe electrocution.

Troubleshooting

If the device does not turn on:

- Make sure that the power cord is plugged properly into the electrical outlet and into the device power socket.

"Smile folds"		5 minutes for each side	Slow up and down motions over the "smile folds"
Chin		5 minutes	Side to side
Jawline		5 minutes for each side	Slow up and down motions along the jaw
Neck		5 minutes for each side	Slow up and down motions across the neck
Decolletage		5 minutes for each 4cm^2 sub-area	Slow circular motions across decolletage
Hands		5 minutes for each 4cm^2 sub-area	Slow circular motions across hands

Recommended Treatment

For long-lasting results, use the device once a week on each desired area for a period of 8 weeks. Maintenance treatments may be performed as needed every 4-8 weeks.

Treatment Area	Recommended Treatment Time	Recommended Device Movement
Forehead 	10 minutes	Slow up and down motion across the forehead
“Eleven” 	5 minutes	Small circular movements or up and down
"Crow's feet" 	5 minutes for each side	Continuous, slow 'side to side' circular motions
Cheeks 	5 minutes for each side	Slow circular motions across the cheek

level button repeatedly until you reach the desired level: 1 or 2.

Note: The device stops the heating when desired temperature is achieved; heating is then re-activated when the temperature beings to drop.

4. Place the treatment fingers on the intended area, ensuring full contact with your skin during the entire treatment. Sensilift will indicate the proper skin contact by the green status light. If Sensilift is not properly placed, the RF energy will not be delivered.

Note: The device will automatically turn-off after 45 seconds of no skin contact.

5. If you do not feel sensations of warmth after approximately 10 seconds, try moving Sensilift on a smaller treatment area, and/or increasing the RF energy.
6. After completing the treatment, wipe off the remaining Base Gel from the treated area.
7. Cleaning instructions for before and after treatment are depicted in the **Cleaning & Storage** section.

Step 3: Post Treatment Care

- We recommend you also use Sensica's anti-aging cream to increase treatment efficacy and ensure long-lasting results (available for purchase at www.sensica.com).
- Apply a high-factor sunscreen (at least 30 SPF) and protect the treated area from sunlight for two weeks following treatment.
- Following the treatment, avoid any usage of hot water on your face and discontinue usage of abrasive or harsh products for two days.
- Makeup may be applied immediately after treatment.

- If the skin in the tested area appears normal and no adverse effects are present, you can proceed with a full treatment at the highest comfortable energy level that you tested.

Note: You may select any RF energy level that you find comfortable. The RF energy level indicates the rate at which the temperature rises when in heating mode. A high RF energy level will enable faster and more even heating. You can achieve maximum results with all RF energy levels, but, for lower levels, it will take more time to reach the effective temperature.

Step 1: Prepare your skin

1. Prior to use, make sure your skin is dry, clean and free of topical products such as: moisturizers, liquids or lint.
2. Remove jewelry from the treatment area as well as from your hands, prior to use.
3. Never use flammable liquids such as alcohol (including perfumes, sanitizers, or other applications containing alcohol) or acetone on the skin one hour before using the device as it may cause burning of the skin.
4. The device is suitable for men. However, male users should make sure their skin is clean-shaven before treatment.

Step 2: Device Operation

Sensilift is equipped with 2 safety sensors which will keep the temperature in the effective and safe range during the entire treatment.

1. Make sure the device is clean, dry and intact prior to use.
2. Apply an even layer of Base Gel, included in your kit, on the area of your skin intended for treatment.
3. Set the Radio Frequency (RF) energy level and the massage speed to the level you selected in your preliminary test treatment. Gently touch the RF energy level selection button ⟆ repeatedly until you reach the desired level: 1, 2, or 3. Gently touch the massage

User Instructions

The device is intended for use on: cheeks, nasolabial "smile" folds, "crow's feet" (next to your eyes), forehead, chin, allowed areas of the neck, décolletage and hands.

What to Expect

During your treatment, it is normal to experience:

- Sensations of warmth created by the RF energy.
- Some mild red or pink color on the skin – during and just after each treatment.

Preliminary Test

The preliminary test is intended to determine your skin's reaction to the treatment and help you find the optimal settings for each treatment area.

Perform a preliminary test on each area you intend to treat 48 hours prior to the first full treatment.

- Prepare your skin as described below in **Step 1: Prepare your skin**.
- Before using the device, apply an even layer of the Base Gel on the treatment area.
- Start your test with energy level 1 and massage level 1 and move the device in small circular motions.
- Increase the energy and massage levels according to your comfort. Stop the test if you do not feel comfortable or if you see adverse effects.
- Perform the preliminary test for a total time of about 3 minutes.
- Wait 48 hours and examine the tested area.


User interface description

1. **ON/OFF and energy selection button** - A gentle touch on the button (after the device is connected to the power supply) will turn the device ON at energy level 1, shown by 1 illuminated light in the energy level display. Once it reaches the highest level, another gentle touch will turn the device OFF. **The device will switch to standby mode when left unattended for 45 seconds.**

2. **Status indicators**

Status indicator	**Blinking yellow light**	The device is in **standby mode**. In case of no light, make sure the power cord is plugged-in properly.
	Steady green light	The device is in **treatment mode**. RF is activated, skin contact is detected, and the treatment is undergoing.
	Blinking alternating green and yellow lights	The device is in **error mode**. The device detected a malfunction and prevents the activation of the RF and massage functions.
RF energy level selection touch button	**1-3 energy levels lights**	When the green light is on, it indicates that the RF energy is active. Touch the RF energy level selection button repeatedly until you reach the desired level: 1, 2, 3 or off.
Massage level selection touch button	**1-2 massage level lights**	Touch the massage level selection button repeatedly until you reach the desired level: 1, 2, or off.

Device Components

Possible Adverse Effects

When using the device according to the Instructions for Use, adverse effects are rare. However, every aesthetic procedure, including the usage of home-use devices, may cause some degree of adverse effects (as listed in the table below).

In case adverse effects appear, stop using the device and consult with your physician if necessary.

Mild adverse effects should be spontaneously relieved within a few hours.

Severity	Frequency	Adverse Effects
Mild	**Common**	Local pain Skin irritation Skin redness Skin sensitivity Allergies or sensitivity to the Base Gel
Moderate	**Rare**	Mild skin burns Pigment changes Excessive redness and swelling Urticaria, purpura or ecchymosis Hematoma
Severe	**Extremely Rare**	Severe burns Infection Bruising Blistering

Any severe incident that has occurred in relation to the device should be reported to us (via customer service).

DRF™ (Dynamic radio frequency) Technology

This device brings a new approach to non-invasive skin rejuvenation featuring Sensica proprietary Dynamic RF technology, combining clinically proven RF (Radio Frequency), commonly used by doctors, with advanced facial massage.

Radio Frequency (RF) is a form of electromagnetic energy, intended to reduce wrinkles, tighten the skin, and renew skin elasticity. It is suitable for all skin tones, easy to use, and provides noticeable, long-lasting results.

RF technology has been used in doctors' clinics for several years and was found to be highly effective. The device has two electrodes that deliver the RF energy into the skin to contract collagen fibers and stimulate the natural production of collagen, which decreases as we age. If used correctly and according to the **Recommended Treatment**, the device will reduce the appearance of wrinkles as well as tighten the skin and renew elasticity. Results may be visible after just one treatment.

Intended Use

The device is intended:

- For women and men who want to: reduce facial wrinkles, tighten their facial skin, and renew skin elasticity.
- For use on: face, cheeks, "smile" folds, "crow's feet", forehead, chin, jawlines, allowed areas of the neck, décolletage, and hands.
- For use on fine to moderate depth wrinkles.
- To be used with the Base Gel. This gel is especially designed to work with Sensica RF devices and ensures proper delivery of RF energy. This will provide the best possible results by improving contact between the device and the skin to improve skin conductivity.

Note: Sensilift is especially designed for your own personal treatment and is intended for a single user and home use only.

genital area, or areas not recommended in this user manual.

- Over the thyroid since it may cause adverse events.
- Over the eyelids, eyes or sides of the eyes (next to your nose), since it may cause eye damage.
- If you have any metal items, such as jewelry or piercings, or tattoos or permanent make up on the designated treatment areas.
- Over sunburns, dark brown or black spots such as: freckles, birth marks, moles or blisters in the treatment area.
- If you are allergic or sensitive to chrome, nickel or any component of the Base Gel.
- Do not use the device over hairy areas of the body, as this may result in the hairs being burnt.

Note: If you have undergone any other aesthetic procedures in the designated treatment area, or if you have treated the skin with a prescribed drug within the past 3 months, or if you have any concerns regarding potential adverse effects, please consult your physician before using the device.

- Have diminished or exaggerated perception of temperature change.
- Have taken Accutane in the past 6 months.
- Have a heat-induced disease or an active skin disease, such as herpes, eczema, rashes etc., in the treatment area.
- Have skin cancer, areas at risk of malignancy, or if you have received radiation therapy or chemotherapy within the past 3 months.

Warnings

Do NOT use the device:

- If you have any permanent implants in the designated treatment area, including but not limited to: silicon or injected chemical materials, metal plates or other metal materials.
- If you had plastic surgery, such as face lifting or eyelid surgery, or if you had aesthetic injections of fillers such as: Hyaluronic Acid, Collagen, fat or Botulinum Toxin, in the treatment area within the past 3 months.
- If you have had a professional chemical peel within the past 3 months with a Hydroxyl Acids (AHA), such as glycolic acid, or β Hydroxyl Acid (BHA), or if you had a professional facial abrasion such as: dermabrasion, facial resurfacing, deep chemical peeling (phenol, Trichloroacetic acid-TCA), and/or ablative laser treatments within the past 3 months.
- Over nostrils, lips of the mouth, ears, breast or

- Before using this device, make sure it is clean, dry and intact. Cleaning instructions are detailed under the Cleaning & Storage section.
- Do not use the device without first applying the Base Gel.
- Do not use the Base Gel after its expiration date.

Note: The Base Gel is for external use only - avoid direct contact with eyes. In the event that such contact does occur, wash the eyes immediately with water.

Contraindications

Do NOT use the device if you:

- Have an active implant, such as a pacemaker or insulin pumps, or if you have any known severe heart disorder, such as arrhythmias.
- Are pregnant or breastfeeding.
- Suffer from severe or uncontrolled endocrine or metabolic conditions, such as diabetes or uncontrolled thyroid disease.
- Have a damaged immune system or autoimmune disorder, such as AIDS or Lupus, or if you use drugs that influence the immune system.
- Suffer from blood clotting disorders, excessive bleeding or bruising, damaged wound healing or if you routinely take blood-thinning medication.
- Suffer from a history of deep vein thrombosis.
- Have a history of collagen disorders or suffer from keloid scarring (growth of extra scar tissue where the skin has healed after an injury).
- Have a general extreme illness, such as hepatitis, liver disease or epilepsy.
- Suffer from sensitive or damaged skin, such as in cases of frequent redness or rash, inflammations, breakouts (such as psoriasis), blisters, sores, infections, acne or extra dry skin.

Carefully read the information below before you begin your first treatment. If you still have any doubts or are unsure that the device is safe for your use, please consult with your physician.

Precautions

- Not to be used by persons under 18 years of age.
- KEEP OUT OF REACH OF CHILDREN.
- The device may be used by persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge provided that they are under supervision or have been instructed in the safe use of the appliance and that they understand the hazards involved.
- Do not place, store or use this device near bathtubs, showers, basins or nearby water, as it may lead to severe electrocution.
- Do not open or repair the device, as it may expose you to dangerous electrical components which may cause serious bodily damage.
- Do not leave the device unattended while plugged into the electricity, as it is a risk of electrocution.
- Do not use the device if any part of the device or its accessories are: not working properly, appear damaged, broken, loose, missing or wet, or if you see or smell smoke.
- Do not use the device with a power supply other than the one provided with the device. This may damage the unit and cause it to stop working. Never use DC adapters which are not approved by the manufacturer.
 In case of power cable damage, the cable should only be replaced by the manufacturer or qualified distributor. Contact our customer service for further assistance.

Table of contents

Precautions 77
Contraindications 76
Warnings 75
Intended Use 73
DRF™ (Dynamic radio frequency) Technology 73
Possible Adverse Effects 72
Device Components 71
User interface description 70
Preliminary Test 69
User Instructions 69
Device Maintenance 64
Device specification parameters 62